

MEK-5DH HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D23YH
	2024-02-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±2.2	25.6	25.5	25.0	26.18
NE%	%	±10.0	65.9	65.9	65.9	66.03
LY%	%	±10.0	15.7	15.7	15.7	17.39
MO%	%	±5.0	10.1	10.1	10.1	7.46
EO%	%	±5.0	7.3	7.3	7.3	8.41
BA%	%	±4.0	1.0	1.0	1.0	0.71
NE	10 ³ /μL	±2.8	16.8	16.7	16.5	17.29
LY	10 ³ /μL	±2.8	4.0	4.0	3.9	4.55
MO	10 ³ /μL	±1.0	2.6	2.6	2.5	1.95
EO	10 ³ /μL	±1.0	1.9	1.9	1.8	2.20
BA	10 ³ /μL	±0.8	0.3	0.3	0.3	0.19
RBC	10 ⁶ /μL	±0.25	5.50	5.46	5.45	5.68
HGB	g/dL	±0.7	17.8	17.5	17.6	17.62
HCT	%	±4.0	55.6	53.7	52.2	62.5
MCV	fL	±8.0	101	98.4	95.8	110.1
MCH	pg	±2.8	32.4	32.1	32.3	31.0
MCHC	g/dL	±3.0	32.0	32.6	33.7	28.2
RDW-CV *	%	±4.0	14.0	14.9	14.2	15.7
PLT	10 ³ /μL	±70	553	544	552	585.3
MPV	fL	±4.0	6.1	7.3	6.2	9.1

Manual Registration	
Lot	23Y
Level	High
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901692M